

PIONEER[®] Density

Automatic Density Measuring Instrument

▶ Automatic Density Measuring Instrument designed for rapid and accurate measurements of density. An American-made product available worldwide delivering good performance at a reasonable cost.



Testing Features



- Small tabletop device
- Analytical balance (available for other measurements)
- Simplified operation

Performance



- Operational Enterprise Software provides data management capability
- Accuracy of $\pm 0.025\%$

Options



N/A (none)

Specifications



Testing Standards:	ASTM D297 Section 16.3 and ISO 2781
Balance Resolution:	0.0001 grams
Density Resolution:	
Weight:	40 lbs (18.1 Kg)
Dimensions:	18" (457mm) L x 10.5" (266mm) W x 30" (762mm) H

